


EXPERIMENT 5
KINETICS ANALYSIS OF AMYLOGUCOSIDASE HYDROLYSIS OF SOLUBLE STARCH


Figure 1: Starch hydrolysis with the final product is glucose (Source: <http://www.sigmaaldrich.com/technical-documents/articles/analytix/enzymatic-food-analysis.html>)

OBJECTIVES:

Student should be able to

1. To understand and explain the enzymatic hydrolysis of starch
2. To explain and calculate enzyme kinetics parameters.

MATERIAL

- Amyloglucosidase enzyme (Enzyme activity >3000 U/mL)
- Dinitrosalicylic acid solution (DNS)
- Sterile distilled water
- Universal bottles (10 each group)
- Glass pipettes
- Pipettor
- Starch (50 g/500 mL = 0.1 g/mL)

METHODOLOGY

- 1) Standard curve for reducing sugar analysis
 By now, each individual group should have their own standard curve. Use this standard curve for today's experiment.

- 2) Enzyme preparation
 Take 1 mL of enzyme solution and dilute it to 25 U/mL using distilled water. This will be the enzyme stock solution.

- 3) Starch solution and enzyme reaction
 - i) Prepare 5 different or varying concentrations of starch solution.
 - a) 1 mL stock solution + 1 mL distilled water
 - b) 1 ml stock solution + 3 mL distilled water
 - c) 1 ml stock solution + 5 mL distilled water
 - d) 1 ml stock solution + 7 mL distilled water
 - e) 1 ml stock solution + 9 mL distilled water
 - ii) Calculate the final concentration of each starch solution.
 - iii) Take 1.0 mL of each concentration and mix with 1.0 mL of enzyme stock solution. Do the same procedures but this time add 1 mL of distilled water instead of enzyme. This will be your control.
 - iv) Incubate the mixture in 55 C waterbath for 15 mins.
 - v) Conduct reducing sugar analysis and estimate the concentration of glucose in each reaction.
 - vi) Plot the Lineweaver-Burk curve and calculate the Michaelis-Menten constant and maximum reaction velocity.

- 4) Different concentration of enzyme
 - i) Add 1.0 mL of starch solution (0.1 g/ mL) in a new bottle and add five different concentration of enzyme in different tubes.
 - a) 1 ml of enzyme + 2ml distilled water
 - b) 1 ml of enzyme + 4 mL distilled water
 - c) 1 ml of enzyme + 6 ml distilled water
 - d) 1ml of enzyme + 8 mL distilled water
 - e) 1 ml of enzyme + 10 mL distilled water
 - ii) Calculate the enzyme concentration of each enzyme solution.
 - iii) Incubate in 55 C for 15 min and estimate the glucose produced.

- iv) Plot an appropriate graph to estimate the glucose concentration for 10 U/mL enzyme activities.

QUESTIONS

1. Explain in details the starch hydrolysis pathway. Include all enzyme involve.